Omnizid

Sprüh- u. Wischdesinfektion lemon, aldehydfrei
Disinfectante spray e disinfectante per superfici profumo al limone, privo di aldeide
Spraydesinfectie en desinfectiedoekjes met citroengreuk, aldehydevrij
Šrodek dezinfekcyjny o zapachu lemon, nie zawierają aldehydów
Rapid Surface Disinfection, aldehyde free
Solution à base d'alcool pour la désinfection rapide des surfaces des dispositifs médicaux, sans aldéhyde
Rozprašovací a stírací desinfekce, bez obsahu aldehydů
Desinfektionsmedel för sprayning och torkning

REF 85322

 OMNIDENT Dental-
Handelsgesellschaft mbH
Gutenbergring 7-9
D-63110 Rodgau
Tel.: +49 (6106) 874-267
BAUA-Nr.: N-29250

LOT

UN 1987

CE 0124

ACHTUNG

e 10 L

100 g contiennent:

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GB

Area of application:

In accordance to 93/42/EEC: Ready-to-use rapid surface disinfectant for non-invasive medical device. In accordance to biocide-directive 98/8/EC: Ready-to-use rapid surface disinfectant for disinfection of smaller surfaces.

PL

Indications for use:

Spray the respective area from a distance of 30 cm until completely covered. Leave to rest. Wipe residue with an uncoloured disposable cloth. Only to be used on alcohol compatible materials (i.e. not suitable for acrylic glass).

Microbiological Effectiveness:

Bactericidal (incl. MRSA) according to VAH/DGHM, EN 1040, EN 1276; levurocidal (Candida albicans) according to VAH/DGHM, EN 1275, EN 1650; tuberculocidal according to VAH/DGHM, EN 14348; membrane covered viruses (i.e. HBV, HIV, HCV) according to BGA(RKI)/DVV, virucidal against Adeno-viruses according to EN 14476.

Effectivity:

Bactericidal (incl. MRSA) in accordance to VAH/DGHM, EN 1276 (high bio load): Conc.: 1 min.

Levurocidal (Candida albicans) in accordance to VAH/DGHM, EN 1650 (high bio load), EN 1275: Conc.: 1 min.

TbB (Tuberculosis)-inactivation in accordance to VAH/DGHM, EN 14348 (high bio load): Conc.: 1 min.

Membrane covered viruses (i.e. HBV, HIV, HCV) in accordance to BGA(RKI)/DVV with bio load: Conc.: 30 sec.

Adeno-viruses (EN 14476, high bio load): Conc.: 30 sec.

100 g contains:

43,0 g Ethanol, 8,0 g 1-Propanol, 0,04 g Didecylidimethylammoniumchlorid.

Safety Instruction:

Flammable liquid and vapour. Keep away from heat/sparks/open flames/hot surfaces. – No smoking. Keep container tightly closed. Wear protective gloves and eye protection. IF ON SKIN (or hair): Remove/take off immediately all contaminated clothing. Rinse skin with water/shower. In case of fire: Use for extinction: CO₂, powder or water spray. Store in a well-ventilated place. Keep container tightly closed. Dispose of contents/container in accordance with local/regional/national/international regulations. Flammable! Do not spray into flame! Flashpoint according to DIN 51755: approx. 25 °C. Max. quantity 50 ml per m². The total quantity per area is not allowed to be more than 100 ml per m² of the ground area. Keep container tightly closed. Packing material should be completely emptied of contents before disposing.

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